



EZSlamX

INSTALLATION/SOFTWARE/USAGE

Q: Question

What are the recommended computer specifications for running EZSlamX software?

A: Answer

For all versions of EZSlamX (3.0.0.0+), we recommend at least the following system requirements:

OS: Microsoft Windows 10 x86/x64 / Windows 11 x64

Processor: Intel Core i5 or equivalent (i7 and newer recommended)

Memory: 8 GB RAM

Storage: 256 GB available space (SSD highly recommended)

USB 2.0 Ports: Minimum 2